

**EUROP-ARM**

Depuis 1973

## Phosphating Armistol 125ml

<https://www.europarm.fr/en/produit-19375-Phosphating-Armistol-125ml>

SKU	Designation	French Law	Vol (ml)	MSRP
EN3282	Phosphating dose for 5L Armistol	Vente libre	125	76.00 € incl. tax

### Concentrated dose of phosphating, to be diluted in 5 liters of water for a tan by boiling.

Product intended for parkerization or the phosphating technique: surface treatment and protection of contemporary and military weapons.

This ARMISTOL product gives weapons a stable appearance and a black-gray color.

The parkerization process is very easy and quick, just one bath is enough.

**Precautions for use:** wearing rubber gloves to avoid touching the parts to be phosphated is recommended.

**Before using the product:** the parts will be white cleaned, without being polished, and carefully degreased. A sandblasted finish is highly recommended.

**Use of the product:** shake the bottle well before pouring its contents into a metal container which will only be used for this purpose (a container made of unenamelled iron and especially not aluminum). Add 5 liters of distilled water (or rainwater).

Provide a small metal box with holes with a handle for small parts and screws.

Carry out phosphating by boiling the parts to be treated for a period which may vary from half an hour to 4 hours.

Check the color (black-gray) by taking out the parts from time to time.

When the desired appearance is obtained, rinse the parts with plenty of water and cover them with vaseline oil, linen or polish.

*Les prix de vente conseillés sont mentionnés à titre indicatif. Les armuriers sont libres de vendre au prix qu'ils souhaitent. Textes et photos non contractuels, sujet à modification.*